

---

# **Timeflux BITalino**

***Release 0.2.1.dev4+g28f7bbb***

**Pierre Clisson**

**Oct 03, 2021**



# CONTENTS

<b>1</b>	<b>Installation</b>	<b>3</b>
1.1	API Reference . . . . .	3
	<b>Python Module Index</b>	<b>7</b>
	<b>Index</b>	<b>9</b>



This plugin provides a driver to connect to [BITalino](#) devices.

When connecting a BITalino device for the first time, the default pairing code is: *1234*.



## INSTALLATION

First, make sure that `Timeflux` is installed.

You can then install this plugin in the *timeflux* environment:

```
$ conda activate timeflux
$ pip install timeflux_bitalino
```

### 1.1 API Reference

This page contains auto-generated API reference documentation.

*timeflux\_bitalino*

---

#### 1.1.1 timeflux\_bitalino

*timeflux\_bitalino.helpers*

---

##### helpers

*timeflux\_bitalino.helpers.transfer*

---

Transfer functions

##### transfer

`timeflux_bitalino.helpers.transfer.VCC = 3.3`

`timeflux_bitalino.helpers.transfer.ECG(signal, resolution=10)`  
ECG value in millivolt [-1.5, 1.5]

`timeflux_bitalino.helpers.transfer.EMG(signal, resolution=10)`  
EMG value in millivolt [-1.64, 1.64]

`timeflux_bitalino.helpers.transfer.EDA(signal, resolution=10)`  
EDA value in microsiemens [0, 25]

`timeflux_bitalino.helpers.transfer.EEG(signal, resolution=10)`  
EEG value in microvolts [-39.49, 39.49]

`timeflux_bitalino.helpers.transfer.PZT(signal, resolution=10)`  
EEG value in percents [-50%, 50%]

`timeflux_bitalino.nodes`

---

## nodes

`timeflux_bitalino.nodes.driver`

---

## driver

`timeflux_bitalino.nodes.driver.TRANSFER`

**class** `timeflux_bitalino.nodes.driver.Bitalino`(*port*, *rate*=1000, *channels*=('A1', 'A2', 'A3', 'A4', 'A5', 'A6'), *sensors*=None)

Bases: `timeflux.core.node.Node`

BITalino driver.

This node connects to a BITalino device and streams data at a provided rate. It is based on the original BITalino Python library, with some performance improvements and careful timestamping.

Two output streams are provided. The default output is the data read from the analog and digital channels. The `o_offsets` output provides continuous offsets between the local time and the estimated device time. This enables drift correction to be performed during post-processing, although no significant drift has been observed during testing.

### Variables

- `o` (*Port*) – BITalino data, provides DataFrame.
- `o_offsets` (*Port*) – Time offsets, provide DataFrame.

### Parameters

- **port** (*string*) – The serial port. e.g. COM3 on Windows; /dev/tty.bitalino-DevB on MacOS; /dev/ttyUSB0 on GNU/Linux.
- **rate** (*int*) – The device rate in Hz. Possible values: 1, 10, 100, 1000. Default: 1000.
- **channels** (*tuple*) – The analog channels to read from. Default: ('A1', 'A2', 'A3', 'A4', 'A5', 'A6').
- **sensors** (*dict*) – The map of attached sensors. If set, transfer functions will be applied. e.g. {"A1": "ECG", "A3": "EMG"}. Default: None.



## Example

```
graphs:

- id: acquisition
  nodes:
  - id: bitalino
    module: timeflux_bitalino.nodes.driver
    class: Bitalino
    params:
      port: /dev/tty.BITalino-02-44-DevB
      rate: 1000
      sensors:
        A1: ECG
        A2: EEG
        A3: EDA
  - id: pub_bitalino
    module: timeflux.nodes.zmq
    class: Pub
    params:
      topic: bitalino
  - id: pub_offsets
    module: timeflux.nodes.zmq
    class: Pub
    params:
      topic: offsets
  edges:
  - source: bitalino
    target: pub_bitalino
  - source: bitalino:offsets
    target: pub_offsets
  rate: 30

- id: display
  nodes:
  - id: subscribe
    module: timeflux.nodes.zmq
    class: Sub
    params:
      topics: [ bitalino, offsets ]
  - id: ui
    module: timeflux_ui.nodes.ui
    class: UI
  - id: debug
    module: timeflux.nodes.debug
    class: Display
  edges:
  - source: subscribe:bitalino
    target: ui:bitalino
  - source: subscribe:bitalino
    target: debug
  rate: 10
```

(continues on next page)

(continued from previous page)

```
# - id: record
#   nodes:
#     - id: sub
#       module: timeflux.nodes.zmq
#       class: Sub
#       params:
#         topics: [ bitalino, offsets ]
#     - id: save
#       module: timeflux.nodes.hdf5
#       class: Save
#   edges:
#     - source: sub:bitalino
#       target: save:bitalino
#     - source: sub:offsets
#       target: save:offsets
#   rate: 1

- id: broker
  nodes:
    - id: broker
      module: timeflux.nodes.zmq
      class: Broker
```

Notes:

**Attention:** Make sure to set your graph rate to an high-enough value, otherwise the device internal buffer may saturate, and data may be lost. A 30Hz graph rate is recommended for a 1000Hz device rate.

Instantiate the node.

**update**(*self*)

Update the input and output ports.

**terminate**(*self*)

Perform cleanup upon termination.

## PYTHON MODULE INDEX

### t

- `timeflux_bitalino`, [3](#)
- `timeflux_bitalino.helpers`, [3](#)
- `timeflux_bitalino.helpers.transfer`, [3](#)
- `timeflux_bitalino.nodes`, [4](#)
- `timeflux_bitalino.nodes.driver`, [4](#)



## INDEX

### B

Bitalino (*class in timeflux\_bitalino.nodes.driver*), 4

### E

ECG() (*in module timeflux\_bitalino.helpers.transfer*), 3

EDA() (*in module timeflux\_bitalino.helpers.transfer*), 3

EEG() (*in module timeflux\_bitalino.helpers.transfer*), 3

EMG() (*in module timeflux\_bitalino.helpers.transfer*), 3

### M

module

timeflux\_bitalino, 3

timeflux\_bitalino.helpers, 3

timeflux\_bitalino.helpers.transfer, 3

timeflux\_bitalino.nodes, 4

timeflux\_bitalino.nodes.driver, 4

### P

PZT() (*in module timeflux\_bitalino.helpers.transfer*), 4

### T

terminate() (*timeflux\_bitalino.nodes.driver.Bitalino*  
*method*), 6

timeflux\_bitalino

module, 3

timeflux\_bitalino.helpers

module, 3

timeflux\_bitalino.helpers.transfer

module, 3

timeflux\_bitalino.nodes

module, 4

timeflux\_bitalino.nodes.driver

module, 4

TRANSFER (*in module timeflux\_bitalino.nodes.driver*), 4

### U

update() (*timeflux\_bitalino.nodes.driver.Bitalino*  
*method*), 6

### V

VCC (*in module timeflux\_bitalino.helpers.transfer*), 3